

Isolating Cells from Construct Protocol

CELLINK Bioink

This is a suggested procedure, please adjust according to your experimental needs.

Protocol aim

The aim of this protocol is to provide instructions for how to isolate cells from a 3D bioprinted and crosslinked construct of CELLINK® Bioink. The cells can subsequently be used for flow cytometry.

Material needed

- Cell-laden 3D bioprinted construct
- PBS with 2 mM EDTA and 0.5% BSA
- 1 ml pipette and pipette tips
- Cell strainer (40 µm nylon)
- 50 mL Falcon tube
- Centrifuge
- Buffer for flow cytometry

1

Protocol

This protocol is for constructs printed in 96-well plate, for other sizes change the recommended volumes accordingly. Constructs of the same condition can be pooled into the same suspension. This protocol is a courtesy of Lisa Oliver, PhD, at the University of Nantes, France.

| Step | Title | Material | Description |
|------|----------------|---|--|
| 1 | Dissociation | <ul style="list-style-type: none">- PBS, EDTA and BSA- Cell laden constructs- 1 ml pipette and pipette tips | <ul style="list-style-type: none">- Prepare the PBS with 2 mM EDTA and 0.5% BSA.- Remove cell culture medium from the construct.- Add 100 µl of buffer to each well.- Gently dissociate the cell-laden construct by pipetting 5-10 times with a 1 ml pipette. <p>Note: Take care not to produce air bubbles .</p> |
| 2 | Cell isolation | <ul style="list-style-type: none">- Cell strainer, 40 µm nylon- 50 mL Falcon tube | <ul style="list-style-type: none">- Place the cell strainer over the centrifuge tube.- Transfer the cell suspension to the cell strainer.- Wash the cell strainer 3 x with 2 ml buffer.- Remove the cell strainer and cap the Falcon tube. |

| | | | |
|---|----------------------------|-----------------------------|--|
| | | - Centrifuge | - Centrifuge the recovered cell suspension at 2000 RPM for 5 minutes. |
| 3 | Flow cytometry preparation | - Buffer for flow cytometry | - Remove the supernatant. - Resuspend the cell pellet in a buffer compatible for flow cytometry and run the flow cytometry analysis according to your protocol of choice. |