

## BIOPRINTING PROTOCOL FOR HUMAN EAR WITHOUT CELLS

---

**Overview:** This protocol is a specific way to a human ear from an MRI scan using CELLINK Start bioink.

**Materials:**

Human ear MRI Scan at a scale of 35% of original size<sup>1</sup>

Slic3r Software (v1.2.9)

CELLINK Start bioink

INKREDIBLE 3D Bioprinter by CELLINK

[Straight tip, 22 GA](#)

[Conical tip, 27 GA](#)

[Conical tip, 25 GA](#)

[Conical tip, 22 GA](#)

**Protocol:**

1. The first step is to upload the human ear MRI to the Slic3r software to create an STL file. Using Slic3r (v1.2.9), convert the 3D model to a bioprinting protocol and toolpath with the following parameters:
  - Layer height = 0.40mm
  - External perimeters extrusion width = 0.45mm
  - Perimeters = 1
  - Infill density = 30%
  - Infill Pattern = Rectilinear
  - Printing speed, F = 600mm/min

Upload the bioprinting protocol with the following name:

*"HumanEar\_Scale35\_LH04\_Infill30\_F600.gcode"*

2. The following bioprinting parameters can be used with the INKREDIBLE 3D Bioprinter by CELLINK using the pneumatic-driven micro-extrusion technology.

- Printing pressure for PH1: 100-110 kPa (Nozzle: [Straight tip, 22 GA](#))
- Printing pressure for PH1: 50-60 kPa (Nozzle: [Conical tip, 27 GA](#))
- Printing pressure for PH1: 30-35 kPa (Nozzle: [Conical tip, 25 GA](#))
- Printing pressure for PH1: 25 kPa (Nozzle: [Conical tip, 22 GA](#))
- Printing speed: 600 mm/min
- Printhead temperature: Room temperature (22°C)
- Printbed temperature: Room temperature (22°C)

### 3. Bioprinting metrics

- a. Time for bioprinting: 5 minutes and 51 seconds per construct
- b. Volume of bioink per construct: 1.2 mL

#### **G-codes:**

*HumanEar\_Scale35\_LH04\_Infill30\_F600.gcode*

#### **Further Information:**

human ear.stl

#### **References:**

N/A