

## FRESH BIOPRINTING PROTOCOL

# PhotoCol<sup>®</sup>

This is a suggested procedure, please adjust it according to your experimental needs. To maintain the sterility of the product, work under sterile conditions.

### Protocol aim

The aim of this protocol is to provide instructions for the bioprinting of complex 3D structures with PhotoCol<sup>®</sup> (Methacrylated Type I Collagen from Advanced BioMatrix) using FRESH printing method. It covers the steps of pre-print procedures, printing, and post-print crosslinking. Changing the parameters in the protocol might change printing conditions such as pressure and speed. This protocol was optimized for the pneumatic Temperature-controlled Printhead installed at the BIO X.

### Materials needed

- PhotoCol<sup>®</sup> (Kit 2,3,4\* – Methacrylated collagen, 20 mM acetic acid, Neutralization Solution, a photoinitiator)
- LifeSupport<sup>™</sup>\*, 1X PBS
- BIO X\*, BIO X6\* 3D Bioprinter
- Cells in suspension and cell culture media
- Positive displacement pipette
- Eppendorf tube with a proper size
- 22G Conical Bioprinting needle (1-inch length)\*
- Well plate or Petri dish
- 3 mL syringes and cartridges with Luer lock connections, female/female Luer lock adaptors\*

\*The product can be purchased in the CELLINK store at [www.cellink.com/store/](http://www.cellink.com/store/).

# Protocol

This protocol is adjusted for printing scaffolds at the final PhotoCol<sup>®</sup> concentration of 6 mg/mL. For other concentrations, recalculations need to be made but the same protocol can be followed. To avoid premature collagen self-assembly, we recommend keeping the collagen and all consumables in the fridge prior to printing.

Once compacted, LifeSupport<sup>™</sup> should be kept in a fridge at 4 °C and used within 12 hours. The ambient temperature should not exceed 23°C during handling or printing.

## 1. Preparation of PhotoCol<sup>®</sup> 10 mg/mL stock solution

### MATERIAL

PhotoCol<sup>®</sup> (lyophilized methacrylated Type I collagen)

20 mM acetic acid

### DESCRIPTION

- In a separate container, weigh the desired mass of PhotoCol<sup>®</sup> for stock solution preparation. Store the remaining PhotoCol<sup>®</sup> material in a fridge.
- Add the corresponding volume of 20 mM acetic acid to the PhotoCol<sup>®</sup> container to achieve the target concentration of 10 mg/mL.

Note. For example, to prepare 3 mL of 10 mg/mL stock solution, add 3 mL of 20 mM acetic acid to 30 mg of lyophilized Photocol<sup>®</sup>.

- Mix on a shaker table or rotator plate at 2-10 °C overnight or until fully solubilized.

## 2. Preparation of LifeSupport<sup>™</sup> bath

### MATERIAL

LifeSupport<sup>™</sup>

1X PBS

Well plate or Petri dish

### DESCRIPTION

- Add 40 mL of cold 1X PBS (4 °C) to LifeSupport<sup>™</sup> tube (sterile powder).

Note: 2 g of sterile LifeSupport<sup>™</sup> corresponding to 15 mL of LifeSupport<sup>™</sup> bath. For detailed directions, visit <https://www.cellink.com/wp-content/uploads/2022/03/FluidformDirectionsforUseVersion6-rev-Sep-2021.pdf>.

- Vortex for 1 minute.
- Put the tube into the fridge (4 °C) for 15 minutes.
- Centrifuge for 5 minutes at 2000 rpm.
- Gently pour off or aspirate the liquid supernatant.
- Grab the tube by the cap, hold it horizontally and gently tap it against a palm 15 times.
- Shake the tube containing dislodged LifeSupport<sup>™</sup> vigorously for 10 seconds. Shake along the length of the tube.
- Centrifuge for further 5 minutes at 2000 rpm.
- The LifeSupport<sup>™</sup> should now be compacted at the bottom of the centrifuge tube. Gently pour off or aspirate any remaining liquid supernatant to leave only the compacted LifeSupport<sup>™</sup> in the bottom of the tube.
- Transfer the resulting LifeSupport<sup>™</sup> bath with a sterile spatula into well plates or Petri dish and store it in a fridge until use.

# 3. Preparation for printing

## MATERIAL

- PhotoCol® 10 mg/mL stock solution
- Neutralization Solution (NS)
- Photoinitiator (PI)
- Ice bath
- Bioprinter (BIO X or BIO X6)
- 3 mL syringes with Luer lock connections
- Cartridge, 3cc
- Luer lock adaptors
- Positive displacement pipette
- Eppendorf tube 1.5 mL
- Cell suspension in cell culture medium of choice
- 22G Conical Bioprinting needle (1-inch length)

## DESCRIPTION

- Place the Temperature-controlled Printhead into the freezer at least 30 minutes before printing.
- Set the BIO X printhead at 5 °C and print bed at 10 °C to guarantee LifeSupport™ bath stability.

Note: Make sure the ambient temperature in the lab is maintained at 21-23 °C, otherwise you may use ice packs inside the printing chamber and on the top of the printer to prevent the printing area from overheating.

- Weigh the necessary amount of PI to achieve a desired concentration of PI in the final bioink.

Note: In our bioinks, we commonly use a PI at 0.25% (w/v) concentration (see Table 1).

- Dissolve the PI in the Neutralization Solution.

Note: Always prepare some extra solution to compensate losses during filtration and transfer from one container to another. For example, dissolve 5 mg of a PI in 90 µL of NS for 1 mL of final bioink.

- Sterile filter the NS/LAP solution using a 0.22 µm filter.
- To prepare 1 mL of 6 mg/mL PhotoCol® solution for printing, transfer 600 µL of PhotoCol® 10 mg/mL stock solution into a sterile Eppendorf using a positive displacement pipette.
- Add 45 µL of NS/LAP solution and use the same pipette to homogenize the resulting solution.
- Add 355 µL of cell suspension with desired cell density and pipette the solution up and down until complete homogenization.

**Table 1.** Suggested values for the preparation of 1 mL of 6 mg/mL PhotoCol® bioink.

$V_{\text{bioink}}$	$C_{\text{final bioink}}$	$C_{\text{stock solution}}$	$V_{\text{stock solution}}$	$V_{\text{NS}}$	$m_{\text{PI}}$	$V_{\text{cell suspension}}$
1 mL	6 mg/mL	10 mg/mL	600 µL	45 µL	2.5 mg	355 µL

- Load a cartridge with the bioink using a pipette.
- Place the cartridge in the printhead and cap with a printing needle.

# 4. Printing

## MATERIAL

- Bioprinter (BIO X or BIO X6)
- Well plate or Petri dish previously filled with LifeSupport™ bath

Cartridge with the PhotoCol® bioink  
22G Conical Bioprinting needle (1-inch length)

#### DESCRIPTION

- Print constructs using suggested parameters:
  - pressure at 8-10 kPa.
  - speed at 3.5 mm/s.

Note: If printability is not as desired, adjust the pressure and/or speed to up/down to extrude more/less material at different speeds.

## 5. Incubation and crosslinking

#### MATERIAL

Cell culture medium

Bioprinter (BIO X or BIO X6) with UV modules of choice for photocuring

#### DESCRIPTION

- Keep the constructs for 10 min at room temperature to ensure initial collagen self-assembly prior to the melting of the supporting bath.
- Incubate the constructs for 30 minutes at 37 °C (5% CO<sub>2</sub> and 95% relative humidity) for further self-assembly of PhotoCol® and LifeSupport™ melting.

Note: Large volumes may require longer times for the supporting bath to fully melt.

- Photocure each construct using the UV module. Time will depend on a chosen PI and its concentration.

Note: It is recommended to use the 405 nm photocuring module for 30 seconds instead of 365 nm if possible, when photocuring PhotoCol® with LAP photoinitiator. Overexposure at the 365 nm wavelength might damage the cells.

- Remove melted LifeSupport™ by replacing it with warm cell media to avoid handling the printed construct. For example, if you printed into a 6-well plate, this can be done by carefully aspirating 2 mL of melted LifeSupport™ out and adding 2 mL of warm cell media. Repeat this process until most of the support bath has been replaced by media.