

## Sectioned Construct

# Immunofluorescence Protocol

Validated for all CELLINK® Bioinks, including the A series, Collagen series, GelMA series, GelX series and CELLINK series. This is a suggested procedure, please adjust according to your experimental needs.

### Protocol aim

The aim of this protocol is to provide instructions for immunofluorescence staining of paraffin embedded, sectioned constructs.

### Material needed

- Microscope slides with sectioned construct according to *Sectioning Protocol*
- Beakers for microscope slides
- Distilled water
- 96% Ethanol
- 100% Ethanol
- Xylene or xylene substitute, e.g. Shandon Xylene Substitute (ThermoFisher, Ref: 9990505)
- PBS/TBS
- Triton 100X or Tween20
- Antigen retrieval buffer (Tri-sodium citrate, distilled water, sodium hydroxide/acetic acid, Triton 100X/Tween20)
- Water bath that can reach 98-100°C
- Super PAP pen (optional)
- Bovine serum albumin (BSA) or serum of host of secondary antibody
- Microscope slide box
- Primary antibody
- Secondary antibody
- DAPI (1 mg/mL)
- Cover glass
- Mounting medium, e.g. Fluoromount G

### Protocol

All handling and use of ethanol and xylene/xylene substitute must be done inside a fume hood with proper PPE and disposed according to local regulation.

Step	Title	Material	Description
1	Prepare antigen retrieval buffer	<ul style="list-style-type: none"> <li>- Tri-sodium citrate</li> <li>- Distilled water</li> <li>- sodium hydroxide/acetic acid</li> <li>- Triton 100X/Tween20</li> </ul>	<p>Note: If antigen retrieval buffer already is prepared, proceed to step 2.</p> <ul style="list-style-type: none"> <li>- Prepare antigen retrieval buffer by:               <ol style="list-style-type: none"> <li>1. Dilute 2.94 g of dihydrate tri-sodium citrate in 1000 mL of distilled water</li> <li>2. Adjust the pH to 6.0 with sodium hydroxide/acetic acid</li> <li>3. Add 0.5 mL of Tween 20/Triton 100x</li> </ol> </li> </ul> <p>Antigen retrieval buffer can be stored 1 month in RT, for extended storage store at 4°C.</p>
2	Prepare blocking solution	<ul style="list-style-type: none"> <li>- BSA or serum of host of secondary antibody</li> <li>- TBS/PBS</li> <li>- Triton 100X/Tween20</li> </ul>	<p>Note: If blocking solution already is prepared, proceed to step 3.</p> <p>Prepare blocking solution by mixing 3% BSA or serum of host of secondary antibody with PBST/TBST (PBS or TBS with 0.05% Triton 100X or Tween20). Store at 4°C.</p>
3	Deparaffination and re-hydration	<ul style="list-style-type: none"> <li>- Microscope slides with sectioned construct</li> <li>- Distilled water</li> <li>- 96% ethanol</li> <li>- 100% ethanol</li> <li>- Xylene or xylene substitute</li> </ul>	<ul style="list-style-type: none"> <li>- Deparaffinize and rehydrate sections by moving microscope slides with sectioned construct through following series:               <ol style="list-style-type: none"> <li>1. Xylene or xylene substitute: 3 x 5 min</li> <li>2. 100% ethanol: 1 min</li> <li>3. 96% ethanol: 1 min</li> <li>4. Distilled water, at least 2 min</li> </ol> </li> </ul>
4	Antigen retrieval	<ul style="list-style-type: none"> <li>- Antigen retrieval buffer</li> <li>- Water bath</li> <li>- Slide container that stays intact at 100°C, e.g. a plastic container</li> </ul>	<ul style="list-style-type: none"> <li>- Put the microscope slides in the container and fill up with antigen retrieval buffer until the slides are covered.</li> <li>- Place the container with samples and retrieval buffer in the water bath.</li> <li>- Turn on heating and let the water reach boiling (98°C-100°C). Let the samples incubate 10 min, starting from when the water starts boiling (above 98°C ok).</li> <li>- Let the samples cool in the buffer.</li> </ul> <p>Note: To speed up cooling put the container with samples and retrieval buffer in cold water after treatment.</p>

5	Preparation of samples	<ul style="list-style-type: none"> <li>- TBS/PBS</li> <li>- Super PAP pen (optional)</li> </ul>	<ul style="list-style-type: none"> <li>- Rinse slides 2 x 3 min in PBS or TBS</li> <li>- Let slides dry and circle the sections with the Super PAP pen to separate/highlight them. This allows for multiple staining on slides.</li> </ul>
6	Blocking	<ul style="list-style-type: none"> <li>- Blocking solution</li> </ul>	<ul style="list-style-type: none"> <li>- Add blocking solution to the sections, incubate 45 min in RT.</li> <li>- Blot of blocking solution, do not rinse.</li> </ul> <p>Note: Always add enough solution to completely cover the sections.</p>
7	Primary antibody	<ul style="list-style-type: none"> <li>- Microscope slide box</li> <li>- Primary antibody</li> <li>- Blocking solution</li> <li>- PBS/TBS</li> </ul>	<ul style="list-style-type: none"> <li>- Make a humified chamber of the microscope box by adding wet paper at the bottom of the box.</li> <li>- Place slides horizontally in the humified chamber.</li> <li>- Dilute primary antibody in blocking solution at recommended concentration, calculate 40-50 <math>\mu</math>L per section.</li> <li>- Add enough primary antibody to completely cover the sections.</li> <li>- Place horizontally in humified chamber overnight at 4°C or in RT 60 min.</li> <li>- Rinse 2 x 3 min with PBS or TBS.</li> </ul>
8	Secondary antibody	<ul style="list-style-type: none"> <li>- Microscope slide box</li> <li>- Secondary antibody</li> <li>- Blocking solution</li> </ul>	<p>Note: <b>Make sure to work without UV in darkness to not bleach the fluorophores.</b></p> <ul style="list-style-type: none"> <li>- Blot off excess PBS/TBS, important to not dilute secondary antibody.</li> <li>- Dilute secondary antibody in blocking solution at recommended concentration, calculate 40-50 <math>\mu</math>l per section.</li> <li>- Add enough secondary antibody to completely cover the sections.</li> <li>- Incubate 60 min in humified chamber at RT.</li> </ul>
9	DAPI stain	<ul style="list-style-type: none"> <li>- DAPI (1 mg/mL)</li> <li>- PBS/TBS</li> </ul>	<p>Note: <b>Make sure to work without UV in darkness to not bleach the fluorophores.</b></p> <ul style="list-style-type: none"> <li>- Dilute DAPI (1 mg/mL) 1:50 in PBS or TBS.</li> <li>- Blot off secondary antibody and add enough DAPI to completely cover the sections.</li> <li>- Incubate 10 min in RT in humified chamber.</li> <li>- Rinse 2 x 3 min with PBS or TBS.</li> </ul>
10	Mount and coverslip	<ul style="list-style-type: none"> <li>- Mounting medium</li> <li>- Cover glass</li> </ul>	<ul style="list-style-type: none"> <li>- Apply a drop of mounting medium to the stained slides.</li> <li>- Cover with a cover glass, apply carefully to avoid air bubbles. Let dry at 4°C horizontally overnight.</li> </ul>