



PEGDA PhotoInk[™] Usage Guide

Purpose

The purpose of this guide is to show first time users how to print with Volumetric's PEGDA PhotoInk[™]. The instructions are written for use with the Lumen X bioprinter and assumes the user has the system and accessories.

Storage

- Recommended storage temperature: 4 °C
- Do NOT freeze the PhotoInk[™].
- Keep in a dry, dark location when not in use.
- Protect from free radical initiators and sources of heat.

Typical Parameters

Power (mW/cm ²)	20	
Layer Height (μm)	100	50
Exposure Time (s)	18	10
1 st Layer Time Scale Factor	4x	4x

Supplies Needed

- 1x Syringe, 1 to 5 mL, depending on the volume required
- 1x Conical nozzle, 22 Ga or 0.400 mm or larger recommended OR for greater precision
- 1x Micropipette, 1000 mL recommended
- 1x Box of tips
- 1x Plastic Razor Blade
- 1x Container, 250 mL or larger, filled with either:
 - o DI water
 - PBS, when printing with cells
- 1x Syringe, optional for clearing channels





• 1x Needle, optional for clearing channels

Before Printing

- 1. Allow PhotoInk[™] to reach room temperature by warming in hand.
- 2. Prepare STL file for printing by progressing through the File, Prepare, and Print tabs.
- 3. Dispense the volume of PhotoInk displayed by LightField. Return the remaining PhotoInk to a dark place close-by if more prints will be conducted or to 4 °C for storage.
- 4. Tap Start.

After Printing

- 1. Remove the build platform from the Lumen X.
- 2. Use the plastic razor blade to gently remove the printed part.
- 3. Place the print in the container of PBS or water to wash the bulk material off.
- 4. If there are channels, a syringe and needle can be used to perfuse the wash solution and remove uncured material.
- 5. If no cells are present, the PBS should be replaced at least three times within 24 hours such that the dye washes away within a day.
- 6. The construct can be sterilized after soaking in ethanol 5 min, followed by soaking in PBS to ensure the ethanol is washed away.