

# Trichrome Staining Protocol

Validated for all CELLINK bioinks, including alginate, nanofibrillated cellulose, collagen and GelMA based bioinks. This is a suggested procedure, please adjust according to your experimental needs.

## Protocol aim

The aim of this protocol is to provide instructions for trichrome staining of sectioned paraffin embedded constructs. Follow *Fixation Protocol*, *Paraffin Embedding Protocol* and *Sectioning Protocol* before starting this protocol.

## Materials needed

- Microscope slides with sectioned construct
- Recommended thickness: 8-12  $\mu\text{m}$
- Beakers for microscope slides
- Distilled water
- 96% Ethanol
- 100% Ethanol
- Xylene or xylene substitute, e.g. Shandon Xylene Substitute (Thermofisher, Ref: 9990505)
- Microscope slide box
- 56°C oven
- Trichrome staining kit, Abcam Cat. No.: ab150686
- Mounting medium, e.g. Vector Laboratories Cat.No.:H-5000
- Cover glass

# Protocol

All handling and use of dyes, ethanol and xylene/xylene substitute must be done inside a fume hood with proper PPE and disposed according to local regulation. If using another Trichrome Staining kit adjust the protocol accordingly.

## 1. Deparaffination and rehydration

### MATERIAL

Microscope slides with sectioned construct

Xylene or xylene substitute

100% Ethanol

96% Ethanol

Distilled water

### DESCRIPTION

- Deparaffinize and rehydrate sections by moving microscope slides with sectioned construct through following series:
  1. Xylene or xylene substitute: 3 x 5 min
  2. 100% ethanol: 1 min
  3. 96% ethanol: 1 min
  4. Distilled water: at least 2 min

## 2. Trichrome stain, Bouin's Fluid

### MATERIAL

Microscope slide box

Bouin's Fluid

Oven at 56°C

### DESCRIPTION

- Make a humidified chamber of the microscope box by adding wet paper at the bottom.
- Blot of the samples and place the microscope slides horizontally in the humidified chamber.
- Add Bouin's Fluid to the sections, close the box and carefully transfer the box to the oven.
- Incubate 60 min followed by a 10 min cooling period.
- Rinse in tap water until sections are completely clear.
- Rinse once in distilled water.

Note: Bouin's Fluid is toxic, keep lid properly closed whenever not within a fume hood and handle with outermost care.

## 3. Trichrome stain, continued

### MATERIAL

Trichrome staining kit

96% Ethanol

## DESCRIPTION

- Apply the staining solutions in following order, always adding enough to completely cover the sections:
  1. Mix equal amounts of Weigert's (A) and Weigert's (B) to maintain working solution of Weigert's Iron Haematoxylin: 5 min.
  2. Rinse in tap water until sections are completely clear.
  3. Rinse once in distilled water.
  4. Biebrich Scarlet/Acid Fuchsin Solution: 15 min.
  5. Rinse in distilled water.
  6. Phosphomolybdic/Phosphotungstic Acid: 15 min.
  7. Remove Phosphomolybdic/Phosphotungstic Acid but do not rinse.
  8. Aniline Blue solution: 7 min.
  9. Rinse in distilled water.
  10. Acetic acid solution (1%): 3-5 min.
  11. Dehydrate slides quickly in two changes of 96% ethanol.

Note: Blot off excess solution between the steps through tapping the slide (on the side, i.e. so it is held vertically) towards a paper covered area of the bench. Especially important after the washes since excess water/ethanol will dilute the staining solutions if not removed.

## 4. Dehydration and clearing

### MATERIAL

100% Ethanol

Xylene or xylene substitute

### DESCRIPTION

- Dehydrate and clear slides by moving the stained microscope slides through following series:
  1. 100% ethanol: 2 x 1 min
  2. Xylene or xylene substitute: 3 x 1 min

## 5. Mount and coverslip

### MATERIAL

Mounting medium

Cover glass

### DESCRIPTION

- Apply a drop of mounting medium to the stained slides.
- Cover with a cover glass, apply carefully to avoid air bubbles.
- Let air dry horizontally overnight.