

## Specification Sheet

# Collagen PREMIUM Neutralization Solution

<b>Product description</b>	The Collagen PREMIUM Neutralization Solution is a solution for pH neutralization of Collagen PREMIUM. For description on how to use with our Collagen PREMIUM product, follow the <b>Neutralization and Casting Protocol Collagen PREMIUM</b> .
<b>Intended use</b>	Solution for pH neutralization of Collagen PREMIUM, <b>Research Grade</b> . For research use ONLY. Not intended for <i>in vitro</i> diagnostics or <i>in vivo</i> uses. Not intended for administration in humans or animals. Produced under sterile and aseptic conditions.
<b>Product number</b>	IKR100000011
<b>Shelf life</b>	12 months, expiration date stated on package.
<b>Storage and handling</b>	Store at 4-8°C. DO NOT FREEZE. Avoid temperature fluctuations.
<b>Safety</b>	Handle in accordance with good hygiene and laboratory safety practices. Read <b>Safety Data Sheet</b> for more information regarding ingredients and potential hazardous compounds.
<b>Related documents</b>	Neutralization Protocol as well as Safety Data Sheet can be downloaded from our website <a href="https://cellink.com/product/collagen-premium">https://cellink.com/product/collagen-premium</a> or scan the QR code below.



Property	Specification	Method
<i>Appearance</i>	Transparent liquid	Visual inspection.
<i>Sterility</i>	Sterile	Tested for the presence of bacteria, fungi and yeast.
<i>Endotoxin level</i>	<0.5 EU/mL	Limulus Amoebocyte Lysate assay, Pharmacopoeia 2.6.14 "Bacterial endotoxins": Method D, accredited by SWEDAC. Accreditation Certification 1240: ISO 15189, 2010-11-22.
<i>pH</i>	7.6±0.2	Assessed with pH meter.